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For more than a quarter of a century, hundreds of thousands of people have used bright light boxes to improve their mood and energy. These light boxes are also called polychromatic or full-spectrum light, because they combine all the colors of visible light to make it appear white. This is evident to anyone who has seen white light pass through a crystal or diamond; the crystal separates the white light into all the colors of the rainbow.

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These polychromatic light boxes (i.e., *many colors*) have been widely accepted as safe and effective. For example, the **Philips Energy Light** devices use these same broad-spectrum fluorescent lamps. To our knowledge, there have been no reported incidents of harm or injury caused by the devices.<sup>1</sup> In fact, long-term studies have shown this type of light use to be safe on the eyes.<sup>2</sup>

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More recently, researchers have discovered that the blue portion of this light is largely responsible for this beneficial reaction. In fact, this color is so effective, researchers refer to it as the “action spectrum of light.”<sup>3</sup> In technical terms, the color is approximately 470 nanometers (nm), and is very much like the color of the summer sky. Studies also show that very low levels of this blue light are as effective as much higher levels of the polychromatic or full spectrum white light.<sup>4 5 6</sup>

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Philips has taken this most effective portion of blue light out of the full-spectrum, white light and has developed the goLITE BLU. As with previous studies, clinical trials with the goLITE have demonstrated it to be as effective as full spectrum white light devices, even though the white light devices are much brighter.<sup>7 8</sup> In comparison to 10,000 lux light boxes, the goLITE BLU only produces 200 lux of light.

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Some people may have heard about the 'blue light hazard', or that blue light may be unsafe for the eyes. This type of blue light is referring to the UV and near UV light, and is on the violet – indigo or deep-blue end of the color spectrum. ***This type of blue light is not the same that is produced by the goLITE BLU.***

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To place this in perspective, it should be remembered that the light from the goLITE BLU is the same type of light found in bright light boxes, and bright light boxes produce 50% – 80% more overall blue than the goLITE. For example, the amount of blue emitted by the goLITE in a study done at Harvard's Brigham and Women Hospital, was approximately 1/5 that of 10,000 lux.<sup>9</sup> As an additional comparison, ambient sunshine produces 95% more blue light than the goLITE BLU.

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To verify the safety of the goLITE device, Philips submitted the goLITE for Ocular Radiation Hazard testing by David H Sliney, Ph.D.<sup>7</sup> The goLITE was determined to be safe with even at a factor of 100 times the normal usage recommendations. In addition, the safety analysis states the following:

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---

\*David H. Sliney, Ph.D. is the current president (2009 – 2010) of the international commission on ocular safety, the Commission International de l'Éclairage (CIE). Dr. Sliney is also a consulting Medical Physicist and head of the Laser/Optical Radiation Program at the US Army Center for Health Promotion and Preventive Medicine.

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“The Philips goLITE meets all of the requirements for an Exempt product (without photobiological risk), since it did not exceed the emission limits provided in both national and international standards for photobiological safety of lamps at the reference test distance of 20 cm.”<sup>10</sup>

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Additionally, in the first published blue LED study, the panels were tested at the lens surface (0 cm distance), and found to be safe:

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
“Although the study anticipated a viewing distance of 50 cm, light safety was assessed at shorter distances as well, including at the panel surface (0 cm). The Food and Drug Administration’s Center for Devices and Radiological Health reviewed the full report and concurred with the analysis and findings, based on the radiological measures provided ....”<sup>11</sup>

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These conclusions are reasonable, considering the amount of blue light emitted from the goLITE is a small fraction (less than 5%) of the blue found in outdoor sunlight. Additionally the goLITE produces no UV or near UV light. Bright white fluorescent lights are considered safe, and the goLITE exposes users to far less total blue light than these 10,000 lux devices.

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**Figure 1:** Spectral distribution comparison of sunlight at 100,000 lux and 20,000 lux (typical summer noon irradiance in direction of our eyes from the horizon) with a traditional 10 000 lux fluorescent bright light (EnergyLight), and with the goLITE. The goLITE delivers substantially less overall blue light than either 10 000 lux fluorescents or sunlight. Additionally, the blue shaded area (representing the spectral distribution of the goLITE Pro) is highly correlated to the 'action spectrum' (460-485 nm) where light is most effective.

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## Quellen

- <sup>1</sup> Partonen, T. Light Therapy, in: Seasonal Affective Disorder; Practice and research. Oxford University Press, 2001. 68-69
- <sup>2</sup> Gallin PF, et al. Ophthalmologic examination of patients with seasonal affective disorder, before and after bright light therapy. Amer J Ophthalmol 1995;119:202-210 (p.208)
- <sup>3</sup> Brainard, G. et al. Action Spectrum for Melatonin Regulation in Humans: Evidence for a Novel Circadian Photoreceptor. Journal of Neuroscience, August 15, 2001, 21(16):6405–6412
- <sup>4</sup> Figueiro, M. et al. Preliminary evidence for spectral opponency in the suppression of melatonin by light in humans. NeuroReport Vol 15 No 2 9 February 2004
- <sup>5</sup> Warman, V et al. Phase advancing human circadian rhythms with short wavelength light. Neurosci Lett 2003; 342:37–40
- <sup>6</sup> Wright, H. et al. Differential effects of light wavelength in phase advancing the melatonin rhythm. J. Pineal Res. 2004; 36:140–144
- <sup>7</sup> Anderson, J. et al. Lux vs. wavelength in light treatment of Seasonal Affective Disorder. Acta Psychiatr Scand 2009;120: 203–212

---

<sup>8</sup> Strong, R. et al. Narrow-Band Blue-Light Treatment of Seasonal Affective Disorder In Adults And The Influence Of Additional Nonseasonal Symptoms. *Depress Anxiety*. 2009;26(3):273-8

<sup>9</sup> Anderson et al. 2009

<sup>10</sup> Sliney, David. Optical Radiation Hazard Analysis; Philips goLITE BLU. 30 August 2009.

<sup>11</sup> Glickman, G. et al. Light Therapy for Seasonal Affective Disorder with Blue Narrow-Band Light-Emitting Diodes (LEDs). *Biol Psychiatry* 2006;59:502–507